SmartCow



an integrated infrastructure for increased research capability and innovation in the European cattle sector

Milk mid-infrared spectra to estimate rumen fermentation parameters

Vanlierde A.¹, Dehareng F.¹, Herremans S.¹, Nichols K.², Kuhla B.³, Eugène M.⁴, Martin C.⁴

¹ CRA-W, Gembloux, Belgium
 ² WUR, Wageningen, The Netherlands
 ³ FBN, Dummerstorf, Germany
 ⁴ INRAE, Theix, France









Rumen parameters




- Interests to measure and monitor rumen parameters
 - Animal health
 - Animal welfare
 - Enconomy
 - Ecology
- Sampling rumen juice and measuring methane emissions
 - Invasive
 - Labourious
 - Expensive analyses or apparatus



Proxies for rumen parameters - Smartcow Focus (WP 6.2)





SmartCow



Development of proxies to estimate rumen parameters → Milk MIR spectra









Proxy for rumen parameters (ammonia, VFA, pH) Milk MIR spectra

Sources of variability in reference datasets

Zootechnical level

- Breed
- Individual variability
- Diet composition
- Diet volume
- **—** ...



Sampling level

- Moment Time before/after feeding
- Sampling location in the rumen?
- Canula/probe/bolus?
- Salivary contamination?
- Standardized apparatus?
- **—** ...






Rumen Juice - Datasets

	CRA-W	INRAE	INRAE	INRAE	FBN
Moment feeding	H0, +2, +4	H-1	H+3	H0, +1.5, +3.5, +5.5, +11.5, +21.5	H-1, +45min
N data	12 * 3	79	28	19 * 6	58 *2
Delay Rumen juice - Milk spectra	D-day	4 days	7-8 days	D-day, -1, +1	D-day, -1, +1
Spectral standardization	Ok	Ok	Ok	To be done	Ok
					1

VFA - pH

Ammonia

Rumen Juice – pH and Volatile Fatty Acids (VFA)



Kaufman, W., H. Hagemeister, and G. Durksen. Adaptation to changes in dietary composition level and frequency of feeding. In: Digestive Physiology and Metabolism in Ruminants, ed. Y. Ruckebusch and P. Thivend. Westport, Ct.: AVI Publishing, 1980, p. 587.




Rumen Juice - pH



- Historical data (INRAE)
 - 3 trials
 - -N = 68
 - Sampling 1 hour **before** eating : N = 44 min : 6.8 ; max : 7.5
 - Sampling 3 hours after eating: N = 24
 min: 6.0; max: 7.8
- Smartcow TNA data (INRAE)
 - 1 trial
 - N = 20
 - Sampling 0, 1.5, 3.5, 5.5, 11.5 hours after feeding : N = 19 for each timing +3.5H \rightarrow min : 6.1; max : 6.6

Poor variability to calibrate No problematic health status to discriminate



Kaufman, W., H. Hagemeister, and G. Durksen. Adaptation to changes in dietary composition level and frequency of feeding. In: Digestive Physiology and Metabolism in Ruminants, ed. Y. Ruckebusch and P. Thivend. Westport, Ct.: AVI Publishing, 1980, p. 587.



Rumen Juice – pH + Milk MIR spectra



Test to develop a proxy based on milk MIR spectra

Sampling 1 hour **before** eating

N	R²c	SEC	R²cv	SECV
44	0.1	0.15	0.0	0.17

No sufficient vairiability in this database

→ To enrich

Sampling 3 hours after eating

N	R²c	SEC	R²cv	SECV
43	0.3	0.33	0.24	0.34

→ Dedicated trials?

→ Standardised protocols!




Rumen Juice - VFA



Consideration of C2, C3, C4, iC4, C5, iC5

 Same conclusion than for pH regarding the variability of reference datasets available

→ Poor variability to calibrate
No problematic health status to discriminate



Kaufman, W., H. Hagemeister, and G. Durksen. Adaptation to changes in dietary composition level and frequency of feeding. In: Digestive Physiology and Metabolism in Ruminants, ed. Y. Ruckebusch and P. Thivend. Westport, Ct.: AVI Publishing, 1980, p. 587.



Rumen Juice - VFA + Milk MIR spectra




INRAE (Historical) mmol/L -1h N = 75

VFA	R²c	SEc	R ² cv	SECV
Acetic	0.34	15	0.16	17
Propionic	0.09	5	0.00	5
Butyric	0.15	3	0.03	3
isoButyrate	0.68	0.2	0.44	0.2
Valeric	0.25	0.9	0.17	0.9
isoValerate	0.68	0.3	0.37	0.5



Rumen Juice - VFA + Milk MIR spectra





$$N = 75$$

 $R^2c = 0.68 - SEc = 0.2$
 $R^2cv = 0.64 - SEcv = 0.2$

N = 75 $R^2c = 0.68 - SEc = 0.3$ $R^2cv = 0.37 - SEcv = 0.5$


- → Not robust but interesting for iC4 and iC5
 - → Consider ratios of VFA
 - → To enrich (health desorders)
 - → Standardised protocols!



Rumen Juice - Ammonia




	CRA-W	INRAE	INRAE	INRAE	FBN
Value	N-NH ₃	N-NH ₃	N-NH ₃	N-NH ₄	NH ₃
Moment feeding	H0, +2, +4	H-1	H+3	H0, +1.5, +3.5, +5.5, +11.5, +21.5	H-1, +45min





Rumen Juice – Ammonia + Milk MIR spectra






- → Interest to consider Milk MIR spectra
- → What is the best timing to consider?
 - → To enrich (health desorders)
 - → Standardised protocols!



Eructed Methane





Eructed Methane + Milk MIR specta

Existing predictive model

CH ₄ ref methods	N data	N cows	SD	R ² cv	SEcv (g/d)
SF ₆ & RC	1,089	299	102	0.68	57

Vanlierde et al. 2020

- Smartcow permitted access to new data (CH₄ (RC) // milk MIR spectra)
 - → Validation step
 - Historical data from FBN
 - n= 233 from 75 dairy cow
 - 392 ± 55 g/d (min: 214 max: 546)
 - TNA data from WUR
 - n = 28 from 28 dairy cows
 - 416 ± 35 g/d (min: 354 max: 475)


- WP5 data from UREAD to be processed
 - 4 dairy cows / 4 periods



Eructed Methane + Milk MIR spectra



Prediction of FBN and WUR data with existing model



$$N = 238$$
, $SD = 51$
 $SEp = 52g CH_4/d$; RPDp = 2
No outliers

→ Error of prediction in line with the model error
 → Detection of high and low emitters



Interest of milk MIR spectra as proxy for rumen parameters but...!



Moments to collect values vary between datasets

Indication of general health status or fermentation profile more than pH or VFA values - no identified « sick » animals in the datasets

- → Too early to conclude especially with datasets available
- → Trial dedicated to this and/or collaborations to enlarge datasets



Standardized protocols!



SmartCow




an integrated infrastructure for increased research capability and innovation in the European cattle sector

Thank you for your attention



SmartCow at a glance





- 11 major RIs distributed in 7 EU countries
- 12 locations, which include 18 installations
- 2500 dairy and 1000 beef cows
- **Networking of RIs** to inventorize resources, harmonize procedures, and share data
- **Joint research activities** to improve experimental methods and phenotyping capability
- **Interaction with stakeholders** to stay in line with industry needs and improve dissemination

http://www.smartcow.eu/stakeholders/

TRAINING PROGRAM

For Scientists, Technicians, Stakeholders, PhD students

- Face-to-face training courses
- Free web-conferences
- One-day study tours in 4 different countries

http://www.smartcow.eu/resources/training/

TRANSNATIONAL ACCESS CALLS

Offers external users (academic and industry) free access to SmartCow RIs

- 30 projects during the 4 years of SmartCow
- Access to around 10,000 cow-weeks

http://www.smartcow.eu/calls/



